


# Anritsu

envision : ensure

ANRITSU INFIVIS CO., LTD.



Solutions for  
Pharmaceutical  
Industries



Weighing Inspection



# Since 1964

Anritsu's inspection systems started with innovation of signal processing.



Checkweigher K501A  
1964



## Anritsu's first Checkweigher was purchased by pharmaceutical company

Anritsu's checkweighers have advanced together with the pharmaceutical industry. We delivered our first checkweigher to a pharmaceutical manufacturer. Learning from the pharmaceutical industry's rigorous attitude toward quality control, Anritsu Infivis has since continued to develop its technologies. Anritsu Infivis is committed to addressing advanced quality assurance issues through closer cooperation with customers and partners in the pharmaceutical industry.

### ▶ 1964




Electronic micrometers

Not so long after desktop electronic calculators were put out into the world, Anritsu gave birth to its checkweigher, which was based on electronic micrometers that electrically measure micron-level displacement using a differential transformer. And that was an innovation. By measuring the difference of weight based on displacement of a coil, we developed a weighcell by ourselves to be built into a machine which is able to reject products - our first checkweigher K501A.



Desktop electric calculator of the time

By MaltaGC at the English language Wikipedia, CC BY-SA 3.0, <https://commons.wikimedia.org/w/index.php?curid=2542893>



Morse printer

### ▶ 1895

In 1895, Anritsu's predecessor Sekisan-sha was founded. The company improved its performance particularly with the production of Morse printers. With a merger in 1908, the company became Kyoritsu Electric Wire, and expanded its operations with the production of public telephones that had just been launched in Japan.



Japan's first public telephone box

Around that time, technologies connected to the current business of the Anritsu Group were discovered; Guglielmo Marconi in Italy succeeded in wireless telegraphy based on radio waves, and Wilhelm Conrad Röntgen in Germany discovered X-rays and succeeded in taking radiograph. Anritsu has a history of advancing together with technologies for signal processing, particularly communications.



[Left] Marconi experimented transatlantic transmission

[Right] Wilhem Röntgen detected X-rays

# Over 80,000 units installed

We disclose six facts in the history of Anritsu for more than half a century.



## Cumulative shipments of checkweigher

► Over 80,000 units as of March 31, 2017

Ever since we shipped our first checkweigher K501A that was put on sale in 1965, we have sold more than 80,000 units of checkweighers in total. We have also sold more than 50,000 units of metal detectors and more than 10,000 units of X-ray inspection systems.



## Countries of destination as of March 31, 2017

► Over 50 countries

Our products are used by customers not only in Japan but also in many countries, from major pharmaceutical manufacturers in the U.S.A., Germany, Italy, China and India, to emerging countries such as Vietnam, Indonesia and Argentina.



## Capsule Checkweigher

► Since 1970

Ever since its early years, Anritsu has developed weight inspection equipment for minute weight products. The K515, capsule checkweigher developed in 1970 featuring five small parallel spring microbalances and Anritsu's proprietary capsule feeding design, had the capacity to weigh 600 capsules per minute at an accuracy of 3 mg.





### Electromagnetic force balance weighcell

▶ Since 1989

Achieved the world's top-level speed of 400 products/minute and accuracy of  $\pm 0.02$  g not by weighing based on gravity but by balancing against electromagnetic force generated by a coil.



### World's first automatic sensitivity setting for metal detectors

▶ Since 1991

Nowadays, most metal detectors are with the function of automatic sensitivity setting, but Anritsu's KD801Ax released in 1991 was the world's first one with the function. With its expertise acquired in the field and total engineering capability in terms of both hardware and software, Anritsu automated phase adjustment that used to be dependent on skilled engineers' intuition and experience.



### In-house development of X-ray inspection system

▶ Since 2000

The KD7203AW is a high-sensitive contaminant detector that generates stable X-rays with an inverter power source and proprietary X-ray tube and uses image processing algorithm developed in-house. Launch of this small-sized X-ray inspection system at an almost half price of conventional large-sized, high-price models facilitated the widespread use of X-ray inspection systems.



# Maximum accuracy: +/-0.5 mg

Anritsu develops all its products, including weighcells, in-house. The results can be seen in the performance of electromagnetic force balance weighcells.



Built-in Multi-Lane  
Weighing System : Weighing unit




## Bringing technologies and wisdoms together

Weighcells, the symbol of Anritsu's originality, become possible only when everything that is required, including hardware, software, production technology, and frontline expertise, is in place. To achieve high speed measurement at the milligram-level accuracy, Anritsu has evolved by making continued efforts for design elaboration and repeating careful simulations, experiments and demonstrations. Anritsu's weighing technologies over half a century is packed into its weighcells.



The weighcell that resistant to vibration noise



Rigid simulation of Weighcell

### ► Dynamic weighing technologies deliver exceptional speed and accuracy

The key to high-speed, high-accuracy weighing is reducing vibration noise and improving response. Weighing accuracy declines with vibration noise generated by non-vertical load applied to the weighcell, which may be caused by dispersion of feeding position of product to be weighed. Anritsu's weighcells have the optimal mechanism for dynamic weighing with improved rigidity due to elaborate analysis of loads in a vertical, horizontal, or torsional direction. Anritsu has also improved response of its weighcells by analyzing weighing vibration to eliminate unnecessary vibration noise because that mixed in weighing signals destabilizes the weight value, which declines processing capacity.

### ► High stability Obtained Class XII, the world's highest standard of EU's Measuring Instrument Directive (MID)

Causes of a decline in weighing stability include a change in room temperature and thermal expansion caused by mechanism elements inside a weighcell. To comply with OIML R51 requirements, Anritsu conducted 3D CAD simulation analysis to halve weight value fluctuations caused by a change in room temperature, compared with its previous model. Anritsu also renewed signal processing to not only to quadruple output resolution compared with its previous model but also to improve rejection accuracy to achieve. Measuring Instruments Directive defines several different classes according to the accuracy of measurement instruments. Anritsu KWS60 series checkweigher is the world's first checkweigher to obtain Class XII\*, the highest standard ever achieved, with its superior weighing properties. Anritsu's new electromagnetic force balance scale has achieved the verification scale interval of  $e=0.05\text{ g}$  (3 to 100 g capacity).

\*as of May 2013, according to our own research



# So many Originals

We disclose six facts that tell Anritsu's technologies.



## Capsule Checkweigher

### ► Stably feed semi-locked capsules

The shuttered magazine, unique shaped retention unit, and ejector allow stable feeding, even for lightweight semi-locked capsules that are difficult to handle, and achieve high-speed, high-accuracy weighing.



## Tablets or Capsules Checkweigher based on sampling

### ► Unique handling technology

Anritsu's tablets or capsules checkweigher suctions tablets or capsules one by one and put them on the weighcell at a certain interval to achieve highly accurate inspection. It is also designed to handle various shapes and sizes of products, such as tablets, capsules,



or spherical Chinese medicine, without replacement of parts.

## Minimum conveyor pitch (Multi-Lane Checkweigher)

### ► 50 mm

Anritsu's small-sized electromagnetic force balance achieves not only the minimum conveyor pitch of 50 mm with proprietary expertise, structural analysis, and stress analysis simulation, but also the maximum accuracy of  $\pm 0.002$  g with improved vibration



resistance and a highly rigid structure.






### X-ray control technology

#### ► Orally disintegrating tablet

Having a structure that allows effective use of low-intensity X-rays, Anritsu's X-ray inspection systems for pharmaceutical products can take high-contrast x-ray images of even pharmaceutical products through which X-ray penetrates easily, making highly sensitive inspection possible.



### Conforming to CFR 21 Part 11

#### ► Only Anritsu in Japan as of June 30, 2018

Anritsu's checkweighers supports authentication, audit trail, and data encryption and decryption, which CFR 21 Part 11 requires. Anritsu allows flexibility in software design because of in-house substrate production.



### Metal sensing technology

#### ► Magnetoreflexion method

Magnetic reflection is a unique detection method based on the principle of metals being magnetized, and invented to inspect metal contaminants inside aluminum packaging materials. Anritsu applied the method to inspect shortage or the number of package inserts printed with magnetic inks.




# Solutions for QC & Productivity

Improves quality and productivity in pharmaceutical manufacturing and R&D.



With increasing population, aging of population, needs for new healthcare, and increasing demand for pharmaceuticals in developing countries, world's demand for pharmaceutical products is increasing at a pace far faster than that of population increase, and expected to increase by about 40% from 2014 to 140 trillion yen (sales value base) in 2020.

Pharmaceutical manufacturers are required to further balance production efficiency and quality maintenance. Anritsu provides a wide variety of high quality solutions to support customers.



Formulation package unit



Filling/dispensing inspection

Tablets



Tablets or Capsules Checkweigher



Capsules



Capsule Checkweigher

Sachets



Built-In Multi-Lane Weighing System

Sheet



Sticks



Multi-Lane Checkweigher



Tube



Dual-Lane Checkweigher



Pouch



Checkweigher for bottles with small diameter



Bottles




Checkweigher for aerosol can



Cans






Topics



Supporting CSV Guidelines: **Validation support**



Anritsu with full validation functions satisfy GMP requirements while supporting accurate and efficient operations. Anritsu also provides IQ/OQ checklists and supports witnessing of PQ to support customers.



Topics



New trends in pharmaceutical inspection: **Non-destructive inside inspection**



**X-ray inspection for 100% inspection of pharmaceuticals:** X-ray inspection systems can conduct shape inspection as to void or broken products, flaw of a tablet, package check, and so on, while inspecting contaminants. The systems can also inspect products in packages that uses aluminum on both sides, and therefore are effective for quality inspection of pharmaceuticals.

**Influence of X-ray radiation on pharmaceuticals:** Jointly with Nagoya City University, Anritsu conducted a survey with regard to influence of X-ray radiation on pharmaceutical quality and confirmed that normal X-ray inspection has no influence on quality.

# Excellent maintenance services

## Global support structure



### Optimize performance at production lines based on the advantage that comes from development and manufacture by ourselves

Anritsu conducts development and manufacture by ourselves, knowledgeable about the characteristics of our products. To make use of those know-how, Anritsu is strengthening our support structure globally.



Our service network covers 50 countries now. Our experienced service engineers provide high-level services.



During service training, we provide know-how on machine settings and trouble shooting in depth so that Anritsu products can perform at the high level that they should be.



Service of sample test being conducted before your purchase of our inspection systems. By using your own product samples, it enables you to confirm exactly how the detection sensitivity can be.



We conduct technical seminars tailored to your requirements. The lecture includes the principle of operation, how to operate and on quality management.

## ANRITSU INFIVIS



Head office

### Head office

#### ANRITSU INFIVIS CO., LTD.

Address : 5-1-1 Onna, Atsugi-shi, Kanagawa-Prf., 243-0032 Japan  
 TEL: +81 46 296 6699 FAX: +81 46 296 6786

Paid-up capital: 1,350 million yen

Sales volume: 17,283 million yen (Year ended March 31, 2017: Non-consolidated)  
 19,588 million yen (Year ended March 31, 2017: Consolidated)

Employees: 404 (as of March 31, 2017: Non-consolidated)  
 533 (as of March 31, 2017: Consolidated)

### Europe, Africa and Middle East

ANRITSU INFIVIS LTD.  
 Luton, United Kingdom

Austria	Luxembourg
Bahrain	Netherlands
Belgium	Norway
Croatia	Oman
Denmark	Poland
Estonia	Portugal
Finland	Qatar
France	Russia
Germany	Saudi Arabia
Greece	South Africa
Hungary	Spain
Ireland	Sweden
Italy	Switzerland
Jordan	Turkey
Kuwait	U.A.E.
Latvia	United Kingdom
Lebanon	Yemen
Lithuania	

### Asia and Oceania

Anritsu Industrial Solutions (Shanghai) Co., Ltd.  
 Shanghai, P.R.China

ANRITSU INFIVIS (THAILAND) CO., LTD.  
 Chonburi, Thailand

Australia	Philippines
China	Singapore
Hong Kong	South Korea
India	Taiwan
Indonesia	Thailand
Malaysia	Vietnam

### Americas

ANRITSU INFIVIS INC.  
 Illinois, U.S.A.

Argentina	Guatemala
Brazil	Honduras
Canada	Mexico
Chile	Panama
Costa Rica	Peru
Ecuador	United States
El Salvador	Venezuela

(as of January 31, 2018)



ANRITSU INFIVIS INC.



ANRITSU INFIVIS LTD.



Anritsu Industrial Solutions (Shanghai) Co., Ltd.



ANRITSU INFIVIS (THAILAND) CO., LTD.



**Anritsu**  
envision : ensure

ANRITSU INFIVIS CO., LTD.

<https://www.anritsu.com/infivis>